



Fast defrosting

Fast defrosting significantly reduces defrosting periods with minimal temperature variations.

MED 720 ULTIMATE

700 I +2 to +15°C USB 🖳 👬 🐼 🗐 🔆

according to DIN

- PRO-ACTIVE-Control: Permanent, proactive monitoring of the performance data and alerts in the case of deviations

 so that you can take countermeasures in good time before a fault threatens your chilled goods; world's most accurate temperature control in refrigerated areas thanks to two standard PT-1000 sensors.
- **Touch-Display,** elegant glass design, capacitive.
- Access control protocol.
- Graphical display of the temperature curve and min/max temperature
- External housing made from galvanised sheet steel (rustproof), with white, anti-scratch powder coating. Length of the plug cable: approx. 3 m.
- Equipped with castors.
- Interior made from smooth aluminium with colourless protective coating. Gibs (every 15 mm) for flexible interior arrangement.
- Interior consists of twelve aluminium drawers with stops, on roller runners. There are two adjustable length dividers and twelve adjustable cross dividers per drawer. Storage surface per drawer = 0.23 m².
- Extra-thick 70-mm insulation, made from high-quality, compression-moulded and environmentally friendly material. Energy saving.

- Self-Closing-Door with easy-to-replace plastic magnetic seal frame, lockable.
- Door hinge on the right-hand side (by default, see illustration), or the left-hand side, can be retrofitted.
- Forced-air cooling, switches off automatically when you open the door, ensures a uniform temperature and minimizes temperature deviations.
- Automatic fast defrosting thanks to time-limited and temperature-monitored reversal of the refrigerant circuit.
- Condensate evaporation in stainless steel heated dish beneath the refrigerator floor.
- Digital data documentation readable via PC-KIT-NET or USB interface and software KIRSCH Datanet.
- Warning functions with visual and audible alarm signal in the case of temperature deviations and power failure. In the case of power failure, the monitoring unit remains in operation for approx. 30 hours on battery power.
- Alarms can be forwarded by network or potential-free contact.
- Antifreeze against sub-zero temperatures.
- Ventilation-enforced refrigerating machine, hermetically sealed, energy saving, low noise, easy to service.

Specification

Capacity	700 litres
Temperature setting	approx. +2 °C to +15 °C
Voltage	220 – 240 V, 50 Hz
Power consumption	250 watts
Normal consumption	1,55 kWh/24 h
Admissible ambient temperature	from +10 °C to +38 °C
Heat emission (max.)	533 watts
Exterior dimensions (in- cluding door handle and distance from wall)	w x d x h = 77 x 98 x 196 cm
Interior dimensions	w x d x h = 62 x 77 x 140 cm (usable depth 10 cm)
Exterior dimensions with door open at 90°	w x d = 77 x 165 cm
Drawer inner dimensions	w x d x h = 57 x 60 x 6 cm
Max. load drawer	24 kg
Weight	Net weight 145 kg, gross weight 180 kg

Optional equipment

- Glass door with lock.
- **LED-Illumination** mounted on side wall.
- Adjustable feet to compense uneven floors.
- Additional length and cross dividers.
- Shelf.
- Condensate container can be emptied manually if no condensation is required (e.g., in an operating theatre).
- 60 Hz cooling machine upon request
- Vario Suport, Storage Shelf on roller runners.
- GSM-Module for alarm text message transmission (e.g., to a mobile phone).
- Temperature documentation:
 - 4-pin PT 100 or PT 1000 temperature sensor, class 1/3 B, including built-in reference body.
 - Mechanical pen-recording thermometer with two waxed paper strips.